



IMPULSE EDITION 400 - CARE BEDS COMPARISON

Features	Impulse Edition 400	Impulse Edition 400 LR
		
<u>Standard Features</u>		
Functions	5 function electric	5 function electric
Backrest and footrest	Standard electric adjustment	Standard electric adjustment
Height adjustment	Standard electric adjustment	Standard electric adjustment
Backrest adjustment	70 °	70 °
Footrest adjustment	30 °	30 °
Height adjustment	25 cm to 82 cm	26 cm to 82 cm
Trendelenburg & reverse	Standard electric adjustment	Standard electric adjustment
Bed dimensions	100 cm (W) x 206 cm (L)	100 cm (W) x 206 cm (L)
Mattress platform	90/100/120 cm (W) x 200 cm (L)	90/100/120 cm (W) x 200 cm (L)
Max user weight	225 kg	225 kg
Fall prevention	Standard	Standard
Transfer position	Standard	Standard
Comfort seating	Standard	Standard
<u>Side rails:</u>		
Split side rail	Aluminium split side rail standard	Aluminium split side rail standard
VGS Split side rail	Aluminium split side rail (one handed) - optional	Aluminium split side rail (one handed) - optional
Full length continues rail	Optional	Optional
Head and footboard	Wooden - designer options available	Wooden - designer options available
Finishing	Various wood finishing options available	Various wood finishing options available
<u>Sleeping surface:</u>		
Metal profile	N/a	N/a
Metal grid	Standard	Standard
ABS	Only available with 90cm width	Only available with 90cm width
Comfort	Optional	Optional
Backrest auto regression	Standard	Standard
<u>Castors:</u>		
Castor type	50 mm twin castors	100/125 mm twin castors
Castor locking	Individual castor locking and axis brake	Central braking
Castors cover	Castors covered	n/a
Auto contour	Standard	Standard
Integrated bed ext.	10/20 cm standard / for full length siderails it is optional	10/20 cm standard / for full length siderails it is optional
Battery backup	Optional	Optional
Guarantee	24 months	24 months
<u>Accessories / Optional</u>		
Lifting/IV pole	Optional	Optional
Bed linen holder	Optional	Optional
Horizontal wall bumper	Optional	Optional
Overbed eating tray	Optional	Optional
CPR emergency lowering	Optional	Optional
Infusion holder for lifting pole	Optional	Optional
Flexible hand control holder	Optional	Optional
LED underbed lighting	Optional	Optional
LED bed lamp	Optional	Optional
Side rail covers	Optional	Optional
Rescue blanket	Optional	Optional