

IMPULSE EDITION 500 - BED SPECIFICATIONS



Features	Impulse Edition 500
<u>Standard Features</u>	
Functions	5 function electric
Backrest, footrest, height	Standard electric adjustment
Backrest adjustment	70 °
Footrest adjustment	40 °
Height adjustment	15 cm to 75 cm
Trendelenburg & reverse	Standard electric adjustment
Bed dimensions	100/120 cm (W) x 205 cm (L)
Mattress platform	90/110 cm (W) x 200 cm (L)
Safe working load	200 kg
Fall prevention	Standard
Transfer position	Standard
Comfort seating	Standard
<u>Side rails:</u>	
Split side rail	Optional
Finishing	Various wood finishing options available
<u>Sleeping surface:</u>	
Metal grid/Metal profile/Comfort	Standard
Backrest auto regression	Standard
<u>Castors:</u>	
Castor type	60 mm twin castors
Castor locking	Central braking
Castors cover	Yes, concealed
Auto contour	Standard
Ultra low position	Standard
Integrated bed ext.	Standard
Battery backup	Optional
Guarantee	24 months
<u>Accessories / Optional</u>	
Lifting / IV pole	Optional
Horizontal wall bumper	Optional
CPR emergency lowering	Optional
LED underbed lighting	Optional
LED bed lamp	Optional
Rescue blanket	N/a