

IMPULSE EDITION XL - BARIATRIC CARE BED



Features	Impulse Edition XL
<u>Standard Features</u>	
Functions	5 function electric
Backrest, footrest, height	Standard electric adjustment
Backrest adjustment	70 °
Footrest adjustment	30 °
Height adjustment	33 cm to 83 cm
Trendelenburg & reverse	Standard electric adjustment
Bed dimensions	110 cm (W) x 206 cm (L)
Mattress platform	100 cm (W) x 200 cm (L)
Max user weight	300 kg
Fall prevention	n/a
Transfer position	n/a
Comfort seating	Standard
<u>Side rails:</u>	
Split side rail	Alluminium split side rail standard
Full length continues rail / VGS	Optional
Finishing	Various wood finishing options available
<u>Sleeping surface:</u>	
Metal grid	Standard
ABS/comfort/metal profile	N/a
Backrest auto regression	Standard
<u>Castors:</u>	
Castor type	150 mm twin castors
Castor locking	Central braking
Castors cover	n/a
Auto contour	Standard
Integrated bed ext.	10/20 cm standard / for full length siderails it is optional
Battery backup	Optional
Guarantee	24 months
<u>Accessories / Optional</u>	
Lifting / IV pole	Optional
Bed linen holder	N/a
Horizontal wall bumper	Optional
Overbed eating tray	Optional
CPR emergency lowering	Optional
LED underbed lighting	Optional
LED bed lamp	Optional
Side rail covers	Optional
Rescue blanket	N/a
Flexible hand control holder	Optional