

Instructions for Use/Assembly Instructions

FS667B/RTS-4”

Raised Toilet Seat with Lid

RAISED TOILET SEAT, DURABLE PLASTIC, 10CM FIXED HEIGHT

Intended use:

The raised toilet seat with a lid is a simple solution for anyone who has trouble getting up from the toilet.

To assemble

- Lift current toilet seat
- Attach Raised Toilet Seat to toilet bowl by securing the seat with clamps on the side of the toilet bowl - Ensure the unit is secured safely and sturdy before use -
- Weight Capacity – 100kg's



White Clamp to secure Toilet Raiser to Toilet Bowl

