

## ***Instructions for Use/Assembly Instructions***

### **SC6040C/SBC**


#### **Swivel Bath Chair with Lock**

##### **Intended use:**

This stainless-steel swivel bath chair is the ideal bath aid, perfect for use by the elderly, or those who suffer from restricted or limited movement.

##### **To assemble:**

- Attach the top part of the Swivel Bath Chair with the Bucket onto the bottom frame by aligning it with the protruding rod on the base
- Use Locking Mechanism by pulling it out gently and releasing it into the grooves on the base and lock gently into desired position
- Weight Capacity – 100kg's



Locking Mechanism to lock Bathchair into desired position

Stabilizer with Lock Nut to secure Bathchair to sides of bath



Top section with bucket seat

Base with protector sleeves

Round black star knob to secure top section onto base – Ensure the top section slides right down into base then tighten knob