

## ***Instructions for Use/Assembly Instructions***

### **SC6040B/SBC**

#### **Swivel Shower Chair with Lock**

##### **Intended use:**

This stainless-steel swivel shower chair is the ideal shower aid, perfect for use by the elderly, or those who suffer from restricted or limited movement.

##### **To assemble:**

- Attach the top part of the Swivel Shower Chair with the Bucket onto the bottom frame by aligning it with the protruding rod on the base
- Release Locking Mechanism by lifting it gently and pressing it down into the hole on the base and lock gently into desired position
- Weight Capacity – 100kg's



Locking Mechanism to allow 360-degree swivel

