



Tel: 021 797 0034 • Fax: 021 797 7664
 Cell & WhatsApp: 076 869 7061
 email: admin@winfar.co.za • website: www.winfar.co.za
 222 Main Road, Wynberg, 7800 - Above Wynberg Pharmacy

Monday - Friday 8am - 5pm
 Closed Saturdays, Sundays & Public Holidays
 Secure Parking in Benjamin Road, behind Wynberg Pharmacy

HOMECARE &
 MOBILITY SOLUTIONS

Please call or email us should you require further
 information on the products in this catalogue.



BASIC STANDARD WHEELCHAIR - 18" BLACK FRAME

CODE	FS809/WCS2.BLK
Footrests	Fixed
Armrests	Fixed
Seat Width	18" - 45 cm
Weight Capacity	100 kg's
Weight of Chair	16,7 kg's

- With full length armrests
- Solid rear wheels
- Black nylon upholstery



BASIC STANDARD WHEELCHAIR - 18" RED FRAME

CODE	FS809/WCS2
Footrests	Fixed
Armrests	Fixed
Seat Width	18" - 45 cm
Weight Capacity	100 kg's
Weight of Chair	16,7 kg's

- With full length armrests
- Solid rear wheels
- Black nylon upholstery



HEAVY-DUTY WHEELCHAIR

CODE	FS975/WCH
Footrests	Detachable
Armrests	Detachable
Seat Width	20" - 50 cm
Weight Capacity	125 kg's
Weight of Chair	18,5 kg's

- With desk length armrests
- Solid rear wheels
- Black nylon upholstery



SUPER HEAVY-DUTY WHEELCHAIR

CODE	FS951B/SUPER-HD
Footrests	Detachable
Armrests	Flip Up
Seat Width	22" - 55 cm
Weight Capacity	150 kg's
Weight of Chair	21,5 kg's

- With padded insert
- Solid rear MAG wheels
- Black nylon upholstery



Wheelchairs Metal Frames

DELUXE WHEELCHAIR WITH QUICK RELEASE WHEELS

CODE	FS908AQ/WCD-QR
Footrests	Detachable
Armrests	Flip Up
Seat Width	18" - 45 cm
Weight Capacity	100 kg's
Weight of Chair	16 kg's

- With drop down back & quick release rear wheels
- Pneumatic rear wheels
- Black nylon upholstery



Wheelchairs - Aluminium Frames

DELUXE LIGHTWEIGHT 18" WHEELCHAIR WITH UNITY BRAKES

CODE	FS908LI/WCL-18"
Footrests	Detachable
Armrests	Flip Up
Seat Width	18" - 45 cm
Weight Capacity	100 kg's
Weight of Chair	13,7 kg's

- With Unity brakes
- Pneumatic rear wheels
- Inclined desk armrests
- Solid Krypton tyres extra optional
- Black nylon upholstery
- Also available in 16" - 40 cm



BEST SELLER

Nappi Code:
282311*001

RECLINING WHEELCHAIR 18"

CODE	FS954GC/RECL
Footrests	Detachable
Armrests	Detachable
Seat Width	18" - 45 cm
Weight Capacity	100 kg's
Weight of Chair	18 kg's

- With high back, elevating footrests & headrest
- Pneumatic rear wheels
- With Unity brakes
- Anti tippers for amputee use
- Solid Krypton tyres extra optional



Nappi Code:
282919*001

ECONOMIC LIGHTWEIGHT WHEELCHAIR

CODE	FS864L/WCL4
Footrests	Fixed
Armrests	Fixed
Seat Width	18" - 45 cm
Weight Capacity	100 kg's
Weight of Chair	12,5 kg's

- With full length armrests
- Pneumatic rear wheels
- Checkered nylon upholstery
- Solid Krypton tyres extra optional



Nappi Code:
282316*001

LIGHTWEIGHT WHEELCHAIR WITH ELEVATING FOOTRESTS

CODE	FS952LCQ/WCL-ELE
Footrests	Elevating & Detachable
Armrests	Flip Up & Adjustable
Seat Width	18" - 45 cm
Weight Capacity	100 kg's
Weight of Chair	17,5 kg's

- Anti tippers for amputee use
- Padded insert
- Solid Krypton tyres extra optional



Nappi Code:
282936*001

Wheelchair Aids

EASY REACHER

- Lightweight aluminium
- Ideal for lifting objects off the floor
- 74 cm long
- User-friendly trigger with magnet
- Trigger works with 1 finger
- Ideal for people with any back / hip issues, restricted reach, mobility issues / a weakened hand grip



(Code: GSOA)

MEMORY FOAM CUSHION

- 55 x 55 cm
- 8 cm high
- Memory foam assist in preventing pressure sores
- Can be used in most chairs & especially wheelchairs
- Removable washable cover



(Code: FS571/MFC-8cm)

Nappi Code:
282921*001

TRANSFER BOARD - STRAIGHT

- Ideal for transferring patients from bed to wheelchair
- Durable lightweight plastic
- 70 cm x 20 cm (Centre thickness 3 cm)



(Code: SC-TWB01/T/Board)

TRANSFER BOARD - CURVED

- Ideal for transferring patients from bed to wheelchair, assisting in getting patients into a car
- The curved shape enables enhanced versatility and user comfort
- Durable lightweight plastic
- Size: 73 cm x 21 cm



(Code: SC-TWB02/T-BOARD)

Wheelchair Repair Service

We offer a Wheelchair, Walker & Rollator Repair Service to our clients



We carry a variety of spares for our Winfar products such as Armrests, Footrests, Front Castors, Rear Wheels, Tyres, Tubes, Backs, Seats, Brake Cables, Handgrips & more.

We also have a Collection / Delivery service in the Cape Metropole
• A collection / delivery charge will apply.

For a FREE quote or enquiries contact Denzil or Lester
Tel: 021 797 0034 or
WhatsApp: 076 869 7061 or
email: admin@winfar.co.za

T's & C's Apply

Mobility Aids

ROLLATOR LOOP BRAKES 8" WHEELS - ALUMINIUM FRAME

BEST SELLER

- Height adjustable handles
- Complete with seat & basket
- Weight capacity 100 kg's (Max)
- Folds up
- Available in Red or Blue

(Code: FS965LH/ROL4)



Napli Code:
282694*001

ROLLATOR WITH 8" WHEELS - ALUMINIUM FRAME & FLIP-UP FOOTRESTS

- Height adjustable handles
- Folds up
- Footrests allow user to be pushed for short distances
- Weight capacity 100 kg's (Max)
- Available in Red or Blue



(Code: FS965LH-B/ROL4)

Napli Code:
1052702*001

TRI-WHEELER RED - BASKET - BAG - TRAY

- Height adjustable handles
- 200 x 50 cm solid castors
- Folds up
- With Bag



(Code: SC5041C/TW-Red)

Napli Code:
1222868*001

KNEE WALKER / SCOOTER

- Height adjustable handles
- Height adjustable & comfortable contoured seat for knee and leg to rest on
- Black Storage Bag
- Rear wheel brakes
- Collapsible handle for easy storage
- Designed to support & assist individuals recovering from foot, ankle, or lower leg injuries



(Code: FS9776/Knee walk)

Napli Code:
1222873*001

FOLDING WALKER WITH 2 WHEELS & PUSH-DOWN BRAKES

- Single Cross bar
- Height adjustable
- Folds up



(Code: SC4010W/FW11)

Napli Code:
1222870*001

FOLDING WALKER (RECIPROCAL)

- Double Cross bar
- Height adjustable
- Folds up



(Code: FS915L/ FWI)

Napli Code:
282890*001

CRUTCH FOREARM BLUE

- A padded cuff fits around the forearm, just above the elbow providing added support and distributing weight evenly
- Preferred for long-term use or when more upper body strength is needed
- Also beneficial for individuals with arthritis
- The grip handle, on the crutch's lower part allows for hand support and control
- The crutch leg is made of aluminium, providing structural support
- Height adjustable and Easy grip handle
- Reflector lights on the handle and back of the cuff
- Anti-slip rubber boot for additional stability and superior grip



(Code: SC1040/CRUTCH FA)

Napli Code:
1222872*001

FOLDING WALKER WITH 4 WHEELS & HAND BRAKES

- 4 castors with cable brakes on rear castors
- Locking brakes on front castors
- Height adjustable
- Folds up



(Code: FS912LJ/WALKER)

Napli Code:
282942*001

QUAD CANE SMALL BASE

- Lightweight aluminium
- Height adjustable
- For left & right hand use
- Offers more stability than a walking stick



(Code: SC3003/QDS)

Napli Code:
1222871*001

QUAD CANE LARGE BASE

- Lightweight aluminium
- Height adjustable
- For left & right hand use
- Offers more stability than the smaller base



(Code: FS921/QDL)

Napli Code:
282700*001

ELBOW CRUTCHES DOUBLE ADJUSTABLE

- Height adjustable handles
- Aluminium lightweight construction
- Double adjustable top to bottom
- Weight capacity 120 kg's



(Code: SC1031/EC-Each)

Napli Code:
282701*001

FOLDING WALKING STICK BRONZE

- Ideal for travelling
- Folds up
- Fits into small travelling bag
- Height adjustable



(Code: SC2022/FWS-BRONZ)

Napli Code:
282702*001

FOLDING WALKING STICK BLACK

- Ideal for travelling
- Folds up
- Fits into small travelling bag
- Height adjustable



(Code: SC2021/FWS-BLACK)

Napli Code:
282702*001

FOLDING WALKING STICK FLORAL DESIGN

- Ideal for travelling
- Folds up
- Fits into small travelling bag
- Height adjustable
- Available in Autumn Brown (SC2141/FWS-FLOR1) Purple & Green (SC2142/FWS-FLOR2)



Napli Code:
282702*001

WALKING STICK T-HANDLE

- Height adjustable
- Aluminium
- Available in Blue Mottled (SC2041/WS-FLORA1) Paisley Blue (SC2041/WS-FLORA2) White & Red (SC2041/WS-FLORA3)



Napli Code:
282702*001

WALKING STICK T-HANDLE

- Height adjustable
- Aluminium
- Available in: Black (SC2031/WS-BLACK) Bronze (SC2032/WS-BRONZE) Silver (SC2033/WS-SILVER)



Napli Code:
282702*001

WALKING STICK CURVE HANDLE

- A comfortable, lightweight mobility aid, with easy-grip curved handle
- Ideal for those in need of just a little extra support when walking
- Made from aluminium
- Available in Black (SC2001/WS BLK CV) Bronze (SC2002/WS BRZ CV)



Napli Code:
282702*001

**WE STOCK A LARGE RANGE
OF WALKING STICKS
& FERRULES**



NB: Designs & Patterns may change from time to time

Commodes

3-IN-1 BASIC COMMODE

BEST SELLER



- Serves as:
Commode,
Raised Toilet Seat
or Safety frame
- Non-folding
- Weight capacity 100 kg's

(Code: FS810/COM10)

Nappi Code:
282435*001

Non-Mobile

4-IN-1 FOLDING COMMODE CHINA IMPORT

- Serves as:
Commode,
Raised Toilet Seat
over toilet or
Safety frame
- Folds up
- Weight capacity
100 kg's



(Code: FS894L/COM1)

Nappi Code:
1010269001

DROP ARM COMMODE CHINA IMPORT

- Non-folding
- Arms drop
for easy
patient transfers
- Height adjustable
- Weight capacity
100 kg's



(Code: FS813/COM/DROP)

Nappi Code:
282437*001

DROP ARM COMMODE TAIWAN IMPORT

- Non-folding
- Arms drop for easy
patient transfers
- Height adjustable
- Weight capacity
120 kg's



(Code: COM7)

Nappi Code:
282437*001

HEAVY DUTY COMMODE

- Weight capacity
200 kg's
- Drop arms
- 23" Seat width
- Height
adjustable



(Code: SC7140HD/COM HD)

Nappi Code:
1222452*001

Mobile

COMMODE FOLDING WITH WHEELS

- Fully foldable commode chair
- It can be used as a bedside commode, raised toilet seat, or safety frame
- Lightweight and easy to move, with locking castor wheels
- Sturdy aluminum frame, Large easy-grip handles
- Height adjustable
- Floor to seat:
Lowest Height: 44 cm
Highest Height: 54 cm
- Hinged plastic toilet seat with 5 litre removable bucket with lid
- Weight Capacity 100 kg's (Max)



(Code: SC7001LW/COMMODE)

Nappi Code:
1222453*001

DROP ARM COMMODE WITH WHEELS

- Lightweight, wheeled commode with drop-down arms for easier transfers
- Built from durable powder-coated metal, it adjusts to fit over toilets, locks safely in place, and includes a removable bucket for simple cleaning
- Stable, Portable, Effortless care
- Height adjustable



(Code: SC7010W/COMMODE)

Nappi Code:
1222454*001

Toilet & Urinary

RAISED TOILET SEAT MOULDED

- Height adjustable
- Suitable for hip-replacement patients
- Weight capacity 100 kg's
- Fits most South African Standard Oval Toilets

Also available with a pu/soft seat,
Code: SC7060K/RTS-PU
Nappi Code: 1222854*001

(Code: SC7060K/RTS)



Nappi Code:
1222864*001

RAISED TOILET SEAT 2-IN-1 S/F/RAME

- Secures over the toilet making sitting down and standing up from the toilet much easier and safer
- Raises toilet seat height by 17 cm
- Designed to fit onto existing toilet bowl, attaches easily & securely to the toilet bowl
- Powder coated metal armrests
- Plastic seat with lid
- Weight capacity: 100 kg's



Nappi Code:
1222865*001

(Code: SC7060H/RTS-SF)

RAISED TOILET SEAT - ARMREST 4" FIXED HEIGHT

- This fixed height, comfortable and secure raised toilet seat, with detachable cushioned aluminum armrests, fits easily atop your existing toilet, increasing safety and ease when using the bathroom.
- Durable plastic seat with lid for hygiene and convenience
- Supports users up to 100 kg's



Nappi Code:
1222866*001

(Code: SC7060B/RTS-4")

RAISED TOILET SEAT 4" WITH LID

- 10 cm (4") fixed height
- Fittings are corrosion free
- Weight capacity 100 kg's
- Fits most South African Standard Oval Toilets



Nappi Code:
282443*001

(Code: FS667B/RTS-4")

BEDPAN SLIPPER WHITE NO LID

- Durable plastic
- Easy to clean



(Code: SC7061B/PBS)

BEDPAN SLIPPER BLUE WITH LID

- Durable blue plastic
- Easy to clean
- Includes lid



(Code: SC7061AC/BPS BLU)

BEDPAN ADULT BLUE

- Durable plastic
- Easy to clean



(Code: SC7061A/BDP BLUE)

MALE URINAL BLUE

- Durable autoclavable plastic lid with 1000 ml capacity



(Code: FS665B/ UN4)

FEMALE URINAL BLUE

- Durable plastic
- Easy to clean
- 900 ml capacity



(Code: FS663B/UN-FEMALE)

URINE BAG STAND

- Designed to hold a urine bag
- Can be placed next to bedside or hang on Hospital bed



(Code: 86258)

MALE URINAL WHITE

- Durable autoclavable plastic lid with 1000 ml capacity



(Code: UN3)

DRAW MAC WITH COTTON SIDES WATERPROOF SOFT PVC

- Cotton tucks allow you to tuck Draw Mac in, to prevent slipping
- Size 90 x 210 cm



(Code: DMCS)

MATTRESS PROTECTOR FITTED WATERPROOF SOFT PVC

- Available in Single, Double, Three Quarter, Queen & King size beds



(Code: MPS)

**We also sell a variety of consumables -
legbags, overnight bags, linen savers,
gloves, d-germ, bandages etc**

Shower Aids

HEIGHT ADJUSTABLE SHOWER CHAIR

BEST SELLER



- Suitable for post hip operation patients
- Handles provide support in raising & lowering from seat
- Weight capacity 100 kg's

(Code: FS7920L/S-CH)

Nappi Code: 282558*001

SHOWER CHAIR WITH BACK BLUE

- Height adjustable
- Detachable foam padded armrests
- Seat with drainage holes
- Weight capacity 100 kg's



(Code: FS7985LA/S-CH)

Nappi Code: 1052705*001

SHOWER CHAIR NO BACK ADJUSTABLE

- Height adjustable
- Non-slip foot grips ensure the shower chair will not slip in the shower
- Armrests to assist with standing
- Removable side arms
- Frame: Aluminum Rust-free
- Seating: Blow moulded plastic
- Weight capacity: 100 kg's



(Code: SC6055AN/S-CH NB)

Nappi Code: 1222446*001

SHOWER BENCH ADJUSTABLE WITH BACK

- Height adjustable
- Moulded seat with drainage holes
- Removable backrest
- Weight capacity 100 kg's



(Code: SC6005Q/SHB-WB)

Nappi Code: 1222447*001

SHOWER BENCH U-SHAPE ADJUSTABLE WITH BACK

- Height adjustable
- U-shaped blow moulded seat for easier bathing
- The two-hole handles in the seat, provide both grip and additional support for the user.
- Weight capacity 100 kg's



(Code: SC6015A/SHB-U)

Nappi Code: 1222448*001

SWIVEL SHOWER CHAIR WITH LOCK

- 360° swivel with locking mechanism
- Height adjustable
- Rust free
- Weight capacity 100 kg's



(Code: SC6040B/SSC)

Nappi Code: 1222449*001

SHOWER CHAIR WALL MOUNTED

- Easy installation
- Mounts & folds up against wall
- Epoxy coated rust proof
- Molded seat with drainage holes
- Two-hole handles in the seat, provide both grip and additional support for the user
- Weight capacity 100 kg's



(Code: SC6190A/SHCW)

Nappi Code: 1222450*001

SHOWER SEAT WALL MOUNTED BRUSHED STAINLESS STEEL

- Brushed stainless steel
- Slatted treated wooden seat
- Includes screws
- Easy to mount
- Weight capacity 100 kg's



(Code: SHCW/F16D)

Nappi Code: 280494*001

SHOWER SEAT WALL MOUNTED BRUSHED WHITE STAINLESS STEEL

- Brushed white stainless steel
- Slatted treated wooden seat
- Includes screws
- Easy to mount
- Weight capacity 100 kg's



(Code: SHCW/F17D)

Nappi Code: 280495*001

HEAVY-DUTY SHOWER CHAIR

- Height adjustable from 38 cm to 55 cm
- Removable backrest
- Ergonomic shaped seat with backrest for extra comfort
- Easy to clean
- Lightweight & easy to transport
- Weight Capacity 180 kg's



(Code: SHOWER/CH-IFIT)

Nappi Code: 186075*001

Bath Aids

SWIVEL BATH CHAIR WITH LOCK

- 360° swivel to assist bathing
- No fitting required
- With positioning lock
- Stainless steel frame
- Fixed armrests for additional user support
- Length from side to side/handle to handle: 70 cm
- Weight capacity 100 kg's



Nappi Code: 282912*001

(Code: SC6040C/SBC)

BATH STOOL ROUND WITH SWIVEL SEAT

- Height adjustable
- NEW Swivel seat top for ease of use
- Suitable for use in bath or shower
- Non-slip tips / ferrules
- Mildew-resistant plastic seat
- Rustproof aluminum frame
- Weight capacity 100 kg's



Nappi Code: 1222451*001

(Code: SC6020D/B-STOOL)

BATH SEAT WITH DRAINAGE HOLES

- Weight capacity 100 kg's
- Drainage holes allows water to drain through seat
- Fits over most baths (width 70 cm)



Nappi Code: 282913*001

(Code: SC6045B/B-SEAT)

BEDSTEP MOULDED PLASTIC SINGLE

- Solid moulded plastic platform
- Aluminium under-frame
- Height adjustable
- Allows easy access to bed - or any other slightly raised surfaces

Length: 41 cm
Width: 23 cm
Height: 15 cm

(Code: SC6050/BSS)

Grab Rails

CHROME GRAB RAILS 12" & 16"

- For extra support & stability in bathroom & toilet
- Non-slip surface in centre provides secure grip
- Available in 30 cm (12") & 40 cm (16")
- Screws included



Nappi Code: 282620*001

(Code: SC6075/GR1-12" & SC6076/GR2-16")

STAINLESS STEEL GRAB RAILS 12" & 16"

- For extra support & stability in bathroom & shower
- Available in 30 cm (12") & 40 cm (16")
- Rust free • Screws included



Nappi Code: 282630*001

(Code: SC6385/GR8-12"-SS & SC6386/GR9-16"-SS)

GRAB RAIL (STUNNING) CRANKED STAINLESS STEEL

- 300 x 300 x 32 mm
- For extra support & stability in bathroom, toilet or shower
- 3 angles allow multiple grips
- Screws included



Nappi Code: 280490*001

(Code: GR/14C)

GRAB RAIL (STUNNING) CRANKED STAINLESS STEEL

- 300 x 300 x 32 mm
- For extra support & stability in bathroom, toilet or shower
- 2 angles for easy gripping
- Screws included



Nappi Code: 280489*001

(Code: GR/13C)

GRAB RAIL FLIP-UP STAINLESS STEEL

- Wall mounted & flips up when not in use
- Toilet roll holder attached
- Screws not included



Nappi Code: 1052715*001

(Code: FS7905S/GR-FLIP)

SAFETY FRAME (STUNNING) STAINLESS STEEL

- 800 x 800 mm
- Wall to floor installation ensures a very secure grip when lowering & raising from the toilet
- Comes in right & left sides
- Screws included



Nappi Code: 1059212*001

(Code: SF/FOTC - RIGHT / SF/FOTC - LEFT)

Models of products may vary from time to time depending on suppliers' availability

Daily Living

Bedroom Aids

COTSIDES

- Fits Single, Double, Three Quarter & Queen size beds



(Code: 2633)

Nappi Code: 525397*001

OVERBED TABLE DUAL TOP LAMINATED WOOD

- Ideal to use when eating, reading or writing
- Swivel top can be used at any angle
- Height adjustable



(Code: FS560/OT)

BACK REST MESH


- Adjusts to 5 positions
- Suitable for beds & leisure use
- Supports back and head
- Soft, padded head cushion included



(Code: SC-W05B/BRST)

BED CRADLE

- Assists in keeping blanket off skin
- Dimensions 55 x 32 x 34 mm
- Essential for bed-bound users and those with sensitive skin or wounds that need to be kept uncovered




(Code: SCBL02/BDCR)

Bedsore Prevention

EGG BOX MATTRESS SINGLE BED

- Helps prevent & relieve bedsore
- High density foam with eggbox shape 1880 mm x 910 mm (Single)
- Also available in Three Quarter, Double, Queen & King size



(Code: PFCS30H00001)

Nappi Code: 1212414*001

MEDICAL SHEEPSKIN

- Helps prevent & relieve bedsore
- Enhances comfort
- Can be handwashed



(Code: SS/XL)

Nappi Code: 115348*001

Sleeping Aids & Support

BED WEDGE WITH COVER

- Developed to provide superior support this wedge cushion provides elevated support for the upper body / legs
- Offers two reclining positions
- Top layer medium density foam with a supportive medium density foam base
- Removable machine washable zip cover, dark blue cotton twill fabric
- Made Locally



(Code: APRODUCTS)

Nappi Code: 121402*001

SLUMBER WEDGE

- Assists in elevating legs or back to aid circulation
- Removable cover



(Code: S/WEDGE)

Nappi Code: 145308*001

Seating

MEDIMAT PARA CUSHION WITH COCCYX CUTOOUT

- High-density foam with eggbox shape
- Most often recommended for wheelchair use as it assists with stability, comfort & pressure relief caused by sitting for long periods
- Made Locally



(Code: PFCS30HPARCOCCYX)

Nappi Code: 1212411*001

MEDIMAT PARA CUSHION PLAIN

- High-density foam with eggbox shape
- Recommended for wheelchair use as it assists with stability, comfort & pressure relief caused by sitting for long periods
- Made Locally



(Code: PFC30HPARA)

Nappi Code: 1212412*001

MEDIMAT WHEELCHAIR EGG BOX CUSHION

- 450 mm x 450 mm
- Suitable for wheelchair & armchairs
- Helps relieve bedsore
- Made Locally



(Code: W/Chair Cushion)

Nappi Code: 1212420*001

SPINE ALIGN WEDGE SEAT

- Cut-out sections relieves pressure on your coccyx & helps reduce back pain
- Made from high density foam for extra comfort & support



(Code: 36967)

Nappi Code: 629428*001

FOAM RING CUSHION COVERED

- Relieves Coccyx pain
- Provides Post-operative & pregnancy pain relief
- High density Latex Foam
- Covered in dark blue stretch trilobal fabric, the cover is removable & machine washable
- Made Locally



(Code: MMRING/LATEXCOV)

Nappi Code: 1212403*001

FOAM RING CUSHION

- Relieves Coccyx pain
- Provides Post-operative pregnancy pain relief
- High Density Latex Foam
- Made Locally



(Code: MMRING/LATEX)

Nappi Code: 1212419*001

HIP REPLACEMENT CUSHION

- Ideal for post hip-operation patients
- High density foam
- Complete with cover & carry strap



(Code: HRP)

Nappi Code: 145297*001

MEMORY FOAM CUSHION CHARCOAL

- High-density memory foam.
- Removable, machine-washable charcoal polycotton cover with zip for easy maintenance
- Made Locally



(Code: MCSPOOOMEM.FOAM)

Nappi Code: 1212410*001

SPINE ALIGN GEL TECH CONTOUR PILLOW

- Contour shape with innovative Gel Tech layer to reduce heat build-up
- Provides for a cool, refreshing sleep experience



(Code: SPINE/GEL-TECH)

Nappi Code: 121402*001

SPINE ALIGN GEL TECH CLASSIC PILLOW

- Contour shape supports your head, neck & upper vertebrae for natural spine alignment
- Reduces tossing & turning whilst sleeping



(Code: SA/GEL-TECH CLASSIC)

Nappi Code: 145308*001

WINFAR MEMORY FOAM CLASSIC PILLOW

- Our Memory Foam Pillows are made from high density visco-elastic memory foam
- Hygienic, hypo-allergenic, anti-bacterial, anti-static & non-toxic
- Made Locally



(Code: PVCLOOOOCLASSIC)

Nappi Code: 1212406*001

WINFAR MEMORY FOAM CONTOUR PILLOW

- Our Memory Foam Pillows are made from high density visco-elastic memory foam
- Hygienic, hypo-allergenic, anti-bacterial, anti-static & non-toxic
- Made Locally



(Code: PVCLOOOOCONTOUR)

Nappi Code: 1212408*001

SPINE ALIGN ERGO BACK SUPPORT

- Moulded to contour the spine
- Provides rigid comfort supporting the spine
- Supports the back & spine muscles from pressure & stress
- Transforms any chair into an ergo chair

(Code: SPINE/ERGO-BACK)

Nappi Code: 1212406*001

POSTURE AID BACK SUPPORT CUSHION

- Ideal for lower back support
- Can be used in car or chair
- Removable machine washable cover, dark blue stretch trilobal fabric
- Made Locally

(Code: PFCS2600000)

Nappi Code: 1212413*001

KNEE RAISER

- Helps elevate knees, aiding blood flow
- Complete with removable machine washable cover, dark blue cotton twill fabric
- Made Locally

(Code: MCSPOOOKNEE.RAIS)

Nappi Code: 1212404*001

SPINE ALIGN MEMORY FOAM CONTOUR LEG PILLOW

- Offers knee support & improves spinal alignment



(Code: SA/CONTOUR LEG)

LUMBAR SUPPORT CUSHION

- Provides comfort and support for the lower back
- Helps ensure better posture
- Soft, yet supportive & flexible
- Made from latex foam
- Covered with a dark blue stretch trilobal fabric
- Made Locally

(Code: MMLUMBLATE)

Nappi Code: 1212405*001

Assistive Devices

BENDABLE SPOON

- Adapted cutlery is perfect for elderly people or those whose condition makes standard cutlery challenging to use
- Made of stainless steel with a built-up black ribbed rubber handle, for improved grip and ease of use
- Dimensions: Length 21 cm, Weight 80 g
- Dishwasher safe

(Code: BENDABLE SPOON)



PLATE GUARD

- Cup on plate guard helps reduce the amount of food spilled when eating allowing an individual to eat independently
- Ideal for plates with a flat edge
- Fits plates 20 cm to 30 cm in diameter
- Top shelf dishwasher and microwave safe
- Easy to detach for cleaning

(Code: PLATEGUARD)



ANEROID BLOOD PRESSURE MONITOR WITH STETHOSCOPE

- Complete manual blood pressure monitoring kit with instructions & travel pouch

(Code: BP2)



FOAM CRUTCH GRIP

- Soft foam grip fits over crutch handle

(Code: FCG)

