



WINFAR is an ISO Compliant Company

For Enquiries: Tel: 087 005 1020 | W/A: 076 869 7061 | admin@winfar.co.za | www.winfar.co.za | 15 Benjamin Road, Wynberg, 7800
Monday - Friday 8am - 4.30pm | Closed Saturday, Sundays & Public Holidays | Secure Parking in Benjamin Street

WHEELCHAIRS: STEEL FRAMES

BASIC STANDARD WHEELCHAIR 18" BLACK FRAME

CODE	FS809/WCS2-BLK
Footrests	Fixed
Armrests	Fixed Full Length
Seat Width	18" / 45 cm
Weight Capacity	100 kg
Weight of Chair	16,7kg



- Solid 8" Front Castors
- Solid 24" Rear Wheels
- Black Powder Coated Frame
- Black Nylon Upholstery

Nappi Code: 282314*001

HEAVY DUTY WHEELCHAIR

CODE	FS975/WCH
Footrests	Detachable
Armrests	Detachable, Desk Length
Seat Width	20" / 50 cm
Weight Capacity	125 kg
Weight of Chair	18,5kg



- Solid 200 x 50 front castors
- Solid 24" rear wheels
- Maroon Powder Coated Frame
- Double Cross Brace
- Black Nylon Upholstery

Nappi Code: 282310*001

BASIC STANDARD WHEELCHAIR 18" RED FRAME

CODE	FS809/WCS2
Footrests	Fixed
Armrests	Fixed Full Length
Seat Width	18" / 45 cm
Weight Capacity	100 kg
Weight of Chair	16,7kg



- Solid 8" front castors
- Solid 24" rear wheels
- Red Powder Coated Frame
- Black nylon upholstery

BEST SELLER

Nappi Code: 282314*001

DELUXE WHEELCHAIR WITH QUICK RELEASE WHEELS

CODE	FS908AQ/WCD-QR
Footrests	Detachable
Armrests	Flip Up, Desk Length
Seat Width	18" / 45 cm
Weight Capacity	100 kg
Weight of Chair	16 kg



- With Drop Back & Quick Release Rear Wheels
- Solid 6" Front Castors.
- Pneumatic 22" Rear Wheels
- Anti Tipppers for Amputee Use
- Double Cross Brace
- Black Nylon Upholstery

Nappi Code: 282930*001

SUPER HEAVY DUTY WHEELCHAIR

CODE	FS951B/SUPER-HD
Footrests	Detachable
Armrests	Flip Up Full Length
Seat Width	22" / 55 cm
Weight Capacity	150 kg
Weight of Chair	21,5 kg



- With padded insert / seat cushion
- Solid 200 x 50 front castors
- Solid 24" rear MAG wheels
- Grey Mottled Powder Coated Frame
- Double Cross Brace
- Anti Tipppers for amputee use
- Black Nylon upholstery

Nappi Code: 282312*001

WHEELCHAIRS: ALUMINIUM FRAMES

RECLINING WHEELCHAIR 18"

CODE	FS954GC/RECL
Footrests	Detachable Elevating
Armrests	Detachable Full Length
Seat Width	18" / 45 cm
Weight Capacity	100 kg
Weight of Chair	18 kg

- With High Back, elevating Footrests & Headrest
- Solid 8" Front Castors
- Pneumatic 24" Rear Wheels
- With Unity Brakes
- Anti Tippers for Amputee Use
- Optional extra: Solid Krypton Tyres

Nappi Code: 282919*001



DELUXE LIGHTWEIGHT 18" WHEELCHAIR WITH UNITY BRAKES

CODE	FS908LI/WCL-18"
Footrests	Detachable
Armrests	Flip Up Desk Length
Seat Width	18" / 45 cm
Weight Capacity	100 kg
Weight of Chair	13,7kg

- With Unity Brakes
- Solid 6" Front Castors
- Pneumatic 22" Rear Wheels
- Solid Krypton Tyres extra optional
- Anti Tippers for Amputee use
- Black nylon upholstery

Nappi Code: 282311*001

BEST SELLER



ECONOMIC LIGHTWEIGHT WHEELCHAIR

CODE	FS864L/WCL4
Footrests	Fixed
Armrests	Fixed Full Length
Seat Width	18" / 45 cm
Weight Capacity	100 kg
Weight of Chair	12,5 kg

- Solid 8" Front Castors
- Pneumatic 24" Rear Wheels
- Checkered Nylon Upholstery
- Optional extra: Solid Krypton Tyres

Nappi Code: 282316*001



EXTRA NARROW LIGHTWEIGHT TRANSIT WHEELCHAIR

CODE	FS804LABJP
Footrests	Flip Up
Armrests	Flip Up
Seat Width	Extra Narrow 14,5" / 37 cm
Weight Capacity	100 kg
Weight of Chair	6,7 kg

- Compact & Portable Lightweight Folding Design
- Lightweight Oxidised Aluminium Frame
- Durable Oxford Nylon Seat & Backrest
- Maintenance-Free Solid Front & Rear Tyres
- Safety Belt Included

Nappi Code: TBC



WHEELCHAIR AIDS

EASY REACHER



- Lightweight aluminium
- Ideal for lifting objects off the floor
- 74 cm long
- 1 Finger, user-friendly trigger with magnet
- Ideal for people with back / hip issues, restricted reach, mobility issues, weakened hand grip

Code: GSOA

WHEELCHAIR TRAY



- Easy clamp-on design for quick attachment & removal
- Durable grey moulded plastic
- Easy to clean
- Lightweight
- Compatible with standard 18 inch wheelchairs

Code: FS505/WCT


WHEELCHAIR SAFETY BELT



- Helps prevent slipping and forward sliding
- Quick release buckle
- Durable woven nylon construction
- Secure screw fit installation
- Compatible with all standard 18 inch wheelchairs

Code: SAFETY BELT

WHEELCHAIR RAMP



- 3 piece modular interlocking ramp kit
- Adjustable height: 20 | 40 | 60 mm
- 750 mm wide ramp surface
- Suitable for indoor and outdoor use
- Lightweight 4 kg aluminium
- Tool-free installation. Quick & easy assembly.
- Weight capacity: 270 kg

Code: WHEELCHAIR/RAMP

WINFAR™ WHEELCHAIR REPAIR SERVICE

We offer a Wheelchair, Walker & Rollator Repair Service for Our Clients

We carry a variety of spares for our WINFAR products, including Armrests, Footrests, Front Castors, Rear Wheels, Tyres, Tubes, Back Seats, Brake Cables, Handgrips & More.

We also offer a *Collection/Delivery Service in the Cape Metropole.

FOR A FREE ASSESSMENT, QUOTE OR ENQUIRIES

Contact Denzil or Lester Tel: 087 005 1020 | W/A 076 869 7061 | E-mail: admin@winfar.co.za

T's & C's Apply | *Nominal Fee Applies



MOBILITY AIDS

ROLLATOR WITH LOOP BRAKES 8" WHEELS | ALUMINIUM FRAME

BEST SELLER



- Height adjustable handles
- Includes Seat and Basket
- Weight capacity: 100 kg (Max)
- Folds Up
- Available in Red or Blue

Code: FS965LH/ROL4

Nappi Code: 282694*001

ROLLATOR WITH FLIP UP FOOTRESTS 8" WHEELS | ALUMINIUM FRAME



- Height adjustable handles
- Folds Up
- Weight capacity 100 kg (Max)
- Footrests allow user to be pushed for short distances
- Available in Red or Blue

Code: FS965LH/B-ROL4BL

Nappi Code: 1052702*001

TRI-WHEELER | BASKET + BAG + TRAY



- Height adjustable handles
- 200 x 50 cm Solid Castors
- Folds Up
- Bag included

Code: 3053

Nappi Code: 282696*001

TRAVEL ROLLATOR WITH FOOTRESTS

- Dual purpose rollator & transport chair
- Foldable steel frame for easy transport & storage
- Comfortable mesh hammock seat
- Height adjustable arm rests
- Flip-up, swing away & detachable footrests
- Storage bag included
- Suitable for indoor / outdoor use
- Max Weight: 100 kg



Code: FS9622/ROL-FRST

ALUMINIUM FOLDING WALKER WITH 4 WHEELS & HAND BRAKES

- 4 Castors with Cable Brakes on Rear Castors
- Locking Brakes on Front Castors
- Height Adjustable
- Folds Up



Code: FS912L J/WALKER

Nappi Code: 282942*001

ALUMINIUM FOLDING WALKER 2 WHEELS & PUSH DOWN BRAKES

- Single Cross Bar
- Height Adjustable
- Folds Up



Code: FS912L/FWI1

Nappi Code: 282697*001

ALUMINIUM FOLDING WALKER RECIPROCAL

- Double Cross Bar
- Height Adjustable
- Folds Up



Code: FS915/FWI

Nappi Code: 282690*001

KNEE WALKER | SCOOTER

- Comfortable alternative to crutches
- Foldable design for convenient transport & storage
- Height adjustable handlebars
- Height adjustable padded knee seat
- Ideal for foot, ankle & lower leg injuries
- Removeable storage basket
- Suitable for indoor & outdoor use
- Weight capacity: 100 kg



Code: 3053

Nappi Code: 257931*001

QUAD CANE | SMALL BASE

- Lightweight aluminium
- Height adjustable
- For left & right hand use
- More stable than a walking stick
- Weight capacity: 100 kg



Code: FS924/QD

Nappi Code: 282699*001

QUAD CANE | LARGE BASE

- Lightweight aluminium
- Height adjustable
- For left & right hand use
- More stable than a walking stick



Code: FS921/QDL

Nappi Code: 282700*001

ELBOW CRUTCHES

- Double Adjustable Top to Bottom
- Height adjustable handles
- Lightweight Aluminium Construction
- Weight Capacity: 120 kg



Code: FS933L/EC

Nappi Code: 282701*001

UNDERARM CRUTCHES

- Height adjustable for individual fit: 112 - 132 cm
- Lightweight Aluminium Construction
- Non-slip rubber ferrules
- Sold in pairs
- Weight Capacity: 100 kg



Code: AUX/CRUTCH-MED

Nappi Code: 186102001

FOLDING WALKING STICK | BRONZE

- Ideal for travelling
- Folds up
- Fits into small travelling bag
- Height adjustable



Code: FS927L/FWS

Nappi Code: 282702*001

FOLDING WALKING STICK | BLACK

- Ideal for travelling
- Folds up
- Fits into small travelling bag
- Height adjustable



Code: FS927L/FWSB

Nappi Code: 282702*001

FOLDING WALKING STICK | FLORAL

- Ideal for travelling
- Folds up
- Fits into small travelling bag
- Height adjustable
- Available in: Autumn Brown: SC2141/FWS-FLOR1 Purple & Green: SC2142/FWS-FLOR2



Nappi Code: 282702*001

FLORAL WALKING STICK | T-HANDLE

- Height adjustable
- Aluminium
- Available in: Blue Mottled: SC2041/WS-FLORA1 Paisley Blue: SC2041/WS-FLORA2 White & Red: SC2041/WS-FLORA3



Nappi Code: 282702*001

WE STOCK A LARGE RANGE OF WALKING STICKS, FERRULES AND FOAM GRIPS.

SHOWER AIDS

HEIGHT ADJUSTABLE SHOWER CHAIR

BEST SELLER

SHOWER CHAIR BLUE WITH BACK & ARMS



- Height adjustable
- Detachable foam padded armrests
- Seat with drainage holes
- Weight capacity: 100 kg

U SHAPE SHOWER CHAIR WITH WHEELS



- U-shaped anti-slip seat for improved hygiene access
- Heavy duty foldable anodized aluminium frame
- Adjustable seat height
- Removeable back rest
- Flip up armrests
- Flip up footrest
- 5-inch castors
- Maximum Weight: 100 kg



- Handles provide support in raising & lowering onto seat
- Suitable for post hip operation patients
- Weight capacity: 100 kg

Code: FS7985LA/S-CH

Nappi Code: 1052705*001

Code: FS962L.SH.FOLD

Nappi Code: TBC

Code: FS7920L/S-CH

Nappi Code: 282558*001

SHOWER SEAT | WALL MOUNTED BRUSHED STAINLESS STEEL

- Brushed stainless steel
- Slatted treated wooden seat
- Includes screws | Easy to Mount
- Weight capacity: 100 kg



SHOWER SEAT | WALL MOUNTED BRUSHED WHITE STAINLESS STEEL

- Brushed white stainless steel
- Slatted treated wooden seat
- Includes screws | Easy to Mount
- Weight capacity: 100 kg



HEAVY DUTY INVACARE SHOWER CHAIR

- Height Adjustable: 38 – 55 cm
- Removeable backrest
- Ergonomic shaped seat with backrest for extra comfort
- Lightweight & easy to transport
- Easy to clean
- Weight capacity: 180 kg



Code: SHCW/F16D

Nappi Code: 280494*001

Code: SHCW/F17D

Nappi Code: 280495*001

Code: SHOWER/CH-IFIT

Nappi Code: 186075*001

BATH AIDS

BATH STOOL ROUND



- Height adjustable 36 - 54 cm
- Suitable for use in bath or shower
- Non-slip tips/ferrules
- Mildew-resistant plastic seat
- Rustproof aluminium frame
- Weight capacity: 100 kg

Code: FS7972L/B-STOOL

Nappi Code: 282916*001

BATH BENCH SEAT/STEP




- Perfect for use as bath step or stool
- Stainless steel frame
- Drainage holes
- Solid & sturdy, with non-slip ferrules
- H: 15cm off the ground | W: 40cm | D: 23cm
- Weight capacity: 100 kg

Code: FS5685/BSS

GRAB RAILS

SAFETY FRAME STAINLESS STEEL



- 800 x 800 mm
- Wall to floor installation ensures a very secure grip when lowering & raising from the toilet
- Comes in right & left sides
- Screws included

Code: SF/FO7C - RIGHT
SF/FO7C - LEFT

Nappi Code: 1059212*001

GRAB RAIL FLIP-UP STAINLESS STEEL



- Wall mounted & flips up when not in use
- Toilet roll holder attached
- Screws not included

Code: FS7905S|GR-FLIP

Nappi Code: 1052715*001

CHROME GRAB RAILS 12" & 16"

- For extra support & stability in bathroom & toilet
- Available in 30 cm (12") & 40 cm (16")
- Screws included



Code: FS7911S-GR/12"
Code: FS7911S-GR-16"

Nappi Code: 282620*001

STAINLESS STEEL GRAB RAILS 12" & 16"

- For extra support & stability in bathroom & shower
- Available in 30 cm (12") & 40 cm (16")
- Rust Free. Screws included



Code: SC6385/GR8-12"-SS
Code: SC6386/GR9-16"-SS

Nappi Code: 282630*001

GRAB RAIL CRANKED STAINLESS STEEL

- 300 x 300 x 300 x 32 mm
- For extra support & stability in bathroom & shower
- 3 angles allow multiple grips
- Screws included



Code: GR/14C

Nappi Code: 280490*001

GRAB RAIL CRANKED STAINLESS STEEL

- 300 x 300 x 32 mm
- For extra support & stability in bathroom, toilet or shower
- 2 angles for easy gripping
- Screws included



Code: GR/13C

Nappi Code: 280489*001

COMMODOES & RAISED TOILET SEATS

3-IN-1 BASIC COMMODE

BEST SELLER

- Serves as Commode, Raised Toilet Seat or Safety Frame
- Removeable bucket with Lid
- Non-Folding
- Weight Capacity: 100 kg



Code: FS810/COM10

Nappi Code: 282435*001

4-IN-1 FOLDING COMMODE

- Serves as Commode, Raised Toilet Seat over toilet or Safety Frame
- Removeable bucket with Lid
- Folds up
- Weight Capacity: 100 kg



Code: FS894/COM1

Nappi Code: 1010269*001

COMMODE FOLDING WITH WHEELS

- Folds up
- Serves as Bedside Commode, Raised Toilet Seat or Safety Frame
- Lightweight. Easy to move
- Locking Castor Wheels
- Sturdy chromed steel frame
- Height adjustable
- 5 Litre removeable bucket with Lid
- Weight capacity 100 kg



Code: FS696/COM1/W

Nappi Code: TBC

RAISED TOILET SEAT MOULDED



- Height adjustable
- Suitable for hip-replacements patients
- Weight capacity: 100 kg
- Fits most South African Standard Oval Toilet

Code: RAISED.T.SEATADJ

Nappi Code: 1222864*001

RAISED TOILET SEAT WITH LID

- 10 cm (4") fixed height
- Fittings are corrosion free
- Weight capacity: 100 kg
- Fits most South African Standard Oval Toilet



Code: FS667B/RTS-4"

Nappi Code: 282443*001

RAISED TOILET SEAT 4" WITH ARMRESTS



- Fixed height, comfortable and secure raised toilet seat, with detachable cushioned aluminium armrests.
- Fits atop your existing toilet.
- Fits toilet with toilet bowl dimensions not less than W*L = 24*30 cm. Max rear lip 4.5 cm.
- Supports users up to 100 kg

Code: SC7060B/RTS-4"

Nappi Code: 1222866*001

URINARY AIDS

MATTRESS PROTECTOR FITTED WATERPROOF SOFT PVC

Sizes: Single, Double, Three Quarter, Queen & King



Code: MPS

DRAW MAC WITH COTTON SIDES WATERPROOF SOFT PVC

- Cotton tucks allow you to tuck Draw Mac in to prevent slipping
- Size: 90 x 210 cm



Code: DMCS

DRIPSTAND MOBILE 5-STAR BASE

- DIMS H:136 – 230 cm
- Top Hooks: includes 4 ram's-horn hooks to hang medical fluid bags safely
- Adjustable Height: Uses a push-pin knob to raise or lower the fluid levels
- Heavy 5-legged base with smooth rolling wheels (castors) preventing tipping while in motion.



Code: FS515/Dripstand

URINE BAG STAND

- Designed to hold a urine bag
- Can be placed next to bedside or hang on Hospital bed



Code: 86258

BEDPLAN SLIPPER BLUE WITH LID

- Durable blue plastic
- Easy to clean
- Includes lid



Code: D115 (01)

BEDPAN ADULT BLUE

- Durable plastic
- Easy to clean



Code: FS661B/BDP

MALE URINAL BLUE

- Durable autoclavable plastic lid
- Capacity: 1000 ml



Code: FS665B/UN4

FEMALE URINAL BLUE

- Durable plastic
- Easy to clean
- Capacity: 900 ml



Code: FS663B/UN-FEMALE

We also sell a variety of consumables -

urine drain bags, leg bags, linen savers, gloves, slide sheets, adult bibs, aprons, d-germ, bandages, needles, syringes, sharps containers etc.

For more information: admin@winfar.co.za

DAILY LIVING

WE ALSO STOCK: A selection of Assistive Devices as well as BPs, STETS etc.

BEDROOM AIDS

COTSIDES




Sizes: Fits Single, Double, Three Quarter & Queen

Code: 2633

Nappi Code: 525397*

EASY GRIP BED RAIL

- Adjustable height bed rail
- Chrome-plated steel frame
- Ultra thin base design
- Extends 3 feet under mattress
- H: 16" W: 20"
- Easy assembly



Code: Bed Rail/Black

Nappi Code: 214451*001

BED CRADLE INVACARE

- Keeps blanket off the skin
- Dimensions: 55 x 32, 5 x 35cm
- Essential for bed-bound users, with sensitive skin or wounds that need to be kept uncovered



Code: 3131

Nappi Code: 143036001

OVERBED TABLE DUAL TOP LAMINATED WOOD



- Ideal to use when eating, reading or writing
- Swivel top can be used at any angle
- Height adjustable

Code: FS560/OT

SEATING SUPPORT

MEDIMAT WHEELCHAIR EGG BOX CUSHION

- 450 mm x 450 mm
- Suitable for wheelchair & armchairs
- Helps relieve bedsores
- Made in SA



Nappi Code: W/Chair Cushion

Nappi Code: 1212420*001

SPINE ALIGN WEDGE SEAT

- Cut-out sections relieves pressure on your coccyx & helps reduce back pain
- Made from high density foam for extra comfort & support



Code: 36967

Nappi Code: 629428*001

MEDIMAT PARA CUSHION WITH COCCYX CUTOUT

- High density foam with eggbox shape
- Recommended for wheelchair users
- Assists with stability, comfort & pressure relief caused by sitting for long periods
- Made in SA



Code: PFC30HP ARCOCCYX

Nappi Code: 12124111*001

MEDIMAT PARA CUSHION PLAIN

- High density foam with eggbox shape
- Most often recommended for wheelchair
- Assists with stability, comfort & pressure relief caused by sitting for long periods
- Made in SA



Code: PFC30HPARA

Nappi Code: 12124112*001

FOAM RING CUSHION COVERED

- Relieves coccyx pain
- Provides post-operative & pregnancy pain relief
- High density Latex Foam
- Covered in dark blue stretch trilobal fabric
- Cover is removable & machine washable
- Made in SA



Code: MMRINGLATEXCOV

Nappi Code: 1212403*001

FOAM RING CUSHION

- Relieves coccyx pain
- Provides post-operative & pregnancy pain relief
- High density Latex Foam
- Made in SA



Code: MMRING/LATEX

Nappi Code: 1212419*001

HIP REPLACEMENT CUSHION

- Ideal for post hip-operation patients
- High density foam
- Complete with cover & carry strap



Code: HRP

Nappi Code: 145297*001

MEMORY FOAM CUSHION CHARCOAL

- High density memory foam
- Removeable & machine washable polycotton cover with zip for easy maintenance
- Made in SA



Code: MCSPO OOMEM.FOAM

Nappi Code: 1212410*001

BEDSORE PREVENTION

EGG BOX MATTRESS SINGLE



- Aids in prevention & relief of bedsores
- High Density Foam, eggbox shape.
- Single: 1880 x 910 mm
- Sizes: King | Queen | Double | 3/4

Code: PFC30H00001

Nappi Code: 1212414*001

MEDICAL SHEEPSKIN



- Aids in prevention & relief of bedsores
- Enhances comfort
- Can be handwashed


Code: SS/XL

Nappi Code: 1193848*001

SLEEPING AIDS & SUPPORT

BED WEDGE WITH COVER

- Elevated support for upper body/legs
- Two Reclining positions
- Top layer medium density foam with a supportive med. density foam base
- Removeable, machine washable zip cover.
- Dark blue cotton twirl fabric
- Made in SA



Code: APRODUCTS

Nappi Code: 121402*001

SLUMBER WEDGE

- Removeable cover
- Assists in elevating legs or back to aid circulation



Code: S/WEDGE

Nappi Code: 145306*001

SPINE ALIGN GEL TECH CLASSIC PILLOW



- Supports your head, neck & upper vertebrae for natural spine alignment
- Reduces tossing & turning whilst sleeping

Code: SA/GEL-TECH CLASSIC

SPINE ALIGN GEL TECH CONTOUR PILLOW



- Contour shape with innovative Gel Tech layer to reduce heat build-up
- Provides a cool, refreshing sleep experience

Code: SPINE/GEL-TECH

SPINE ALIGN ERGO BACK SUPPORT

- Moulded to contour the spine
- Provides rigid spine support and comfort
- Supports the back & spine muscles from pressure & stress
- Transforms any chair into an ergo chair



Code: SPINE/ERGO-BACK

POSTURE AID BACK SUPPORT CUSHION

- Ideal for lower back support
- Can be used in car or chair
- Removeable, machine washable, dark blue, stretch trilobal fabric cover
- Made in SA



Code: PFC2600000

Nappi Code: 121413*001

WINFAR MEMORY FOAM CLASSIC PILLOW



- Made from high density visco-elastic memory foam
- Hygienic, hypo-allergenic, anti-bacterial, anti-static, non-toxic
- Made in SA

Code: PVCL0000CLASSIC

Nappi Code: 1212406*001

WINFAR MEMORY FOAM CONTOUR PILLOW



- Made from high density visco-elastic memory foam
- Hygienic, hypo-allergenic, anti-bacterial, anti-static, non-toxic
- Made in SA

Code: PVCL0000CONTOUR

Nappi Code: 1212408*001

LUMBAR SUPPORT CUSHION

- Provides comfort & support for lower back
- Helps ensure better posture
- Soft yet supportive & flexible
- Made from Latex Foam
- Cover: Dark blue stretch trilobal fabric
- Made in SA



Code: MMLUMBLATE

Nappi Code: 121405*001

KNEE RAISER

- Elevates knees aiding blood flow
- Cover: Dark Blue, Twirl fabric. Removeable & machine-washable.
- Made in SA



Code: MCSPOO OKNEE.RAIS

Nappi Code: 1212404*001

SPINE ALIGN MEMORY FOAM CONTOUR LEG PILLOW



- Offers knee support
- Improves spinal alignment

Code: SA/CONTOUR LEG